

SimpleMed Suit

Special medication dispenser

A customized compliance solution for any medication modality. Designed to combine patient's needs with practicality & ease of use; thus making it the most effective solution for improving medication adherence with special meds. SimpleMed Suit makes use of auditory & visual alerts to remind a patient to take medication at the designated time.

Optimizes service provider - caregiver - patient connection by using Vaica's cloud.



How does SimpleMed Suit work?



Once activated, SimpleMed Suit automatically reads its settings from Vaica's cloud.



Slightly before the scheduled medication time, the LED lights up.



Service provider pre-programs specific device settings via Vaica's portal.



When the scheduled time arrives an audible alert is sounded.



7.

6.

If the patient takes the med on time, a message will be sent to Vaica's cloud.

If the patient doesn't take the med on time, an e-mail / SMS will be sent automatically to the patient, service providers and /or caregivers.

Upon each extraction of medication, the unique medication ID will be sent to Vaica's cloud.

SimpleMed Suit Features:

- GSM connectivity cloud services based on an internal cellular modem
- Additional reminders can be programmed (ex: check blood sugar / apply eye drops)
- Adjustable schedule Personalized to fit patient's needs
- Remote Programming Configuration can be done through a secure web access
- Built-in hub* for vital sign monitoring devices (*optional)
- Built-in NFC reader* enables unique med ID identification (*optional)

Dimensions: Customizable | Temperature: -10°C to 45°C (-50°F to 110°F)

Main Power: 120-240 VAC, 50/60 Hz | Backup Power: Li-ion Battery, Up to 10 day backup Communication: EGSM 850/900/1800/1900 MHz | Bluetooth: Bluetooth 4.1* optional

